



A regenerable VOC control system (RVCS) for characterizing properties of sorbents used in separation technologies

Sara D. Nolek

Oscar A. Monje

ESC - Team QNA

Space Life Sciences Laboratory

Kennedy Space Center, FL 32899

Team
QinetiQ
North America

© QinetiQ North America, Inc. 2010

Introduction



Bed and System Design

Methods

- Thermal Regeneration
- Static Adsorption Test
- Pressure Swing Adsorption and Regeneration

Testing


Discussion



- Adsorptive Capacities
- Roll-over Effect

Conclusions

Bed design


- 2 bed diameters
 - Adjust contact times with gas stream flowrate and bed diameter
- Measurements
 - Temperature
 - 25%, 50%, 75% length of bed
 - Inlet/outlet
 - Humidity
 - Inlet/outlet
 - Pressure
 - Inlet/outlet







System Design

- Gas stream
- Operation through test bed
- Analysis
 - FTIR
 - GC
 - Dewpoint hygrometer
 - CO₂, P, T, Rh sensors
 - Datalogger





Gas Stream






Analysis





Methods cont.

- Static adsorption test
 - Simulated spacecraft gas stream
 - 6000 ppm CO₂
 - 52% Relative Humidity
 - 1-5 ppm VOC
 - 23°C ambient temperature




Typical data set collected in the RVCS during a static test measuring breakthrough curves of CO₂ and water vapor

Methods


■ Thermal regeneration

- Zeolites require high temperature for desorption ($>300^{\circ}\text{C}$)
- Heating system designed/built for high temperature but also ability to hold stable at intermediate temperature
- Reverse flow sweep gas



Methods cont.

- PSA (Pressure Swing Adsorption)
 - <0.5 Torr (0.04 kPa) vacuum
 - Capability to cycle at integral minute intervals



Inlet (PreRH) and outlet (PostRH) relative humidities of a Zeolite 13X bed undergoing 8 min PSA cycles. The bed kept the air exiting the bed dry.

Testing



- Flow rate/carrier gas/pre & post accuracy
 - Does N₂ compete with the CO₂ or the VOC adsorption?
 - Test with He vs. N₂ as carrier gas
 - Does the flow rate affect the CO₂ breakthrough or VOC adsorption?
 - Test with higher flow vs. lower flow
 - Is the carryover from previous measurement enough to skew results?
 - Test “stickiness” of VOC for pre/post concentration error

Testing cont.




- What are the effects of
 - Single VOC on CO₂ adsorption/CO₂ on VOC adsorption?
 - Water on CO₂ adsorption?
 - Water on VOC adsorption/VOC on water adsorption?
 - Adsorption effects of multiple VOCs?
- Comparison with adsorption models

Discussion



- Adsorptive capacities
- Rollover effect
- Stickiness of VOCs to system components
- Small diameter sorbents (powder)

Adsorptive capacity and rollover effect



Breakthrough curve of toluene + water. Once water begins to breakthrough, roll-over of toluene is observed as the water displaces the adsorbed toluene off the bed.

Conclusions



- Recap of system uses
- Correlation with models
- Future capabilities to be added

Acknowledgement



- Support of this research was provided by NASA's Exploration Life Support and Habitation Systems Project– Atmosphere Revitalization Element and NASA's ISS research.